

Transfection into Mouse-Rat: CULTURED EMBRYOS by Ex Utero Electroporation

The NEPA21 is the only device on the market to approach Electroporation from the perspective of optimising delivered energy.

- Compared to devices from other suppliers, the NEPA21 system offers the researcher a level of previously unavailable control over energy delivery to the electroporation target. This control is generated via unique electroporation pulse-output configurations, client-confirmed protocols and application-customised electrodes.
- With this market-leading control and (user-independent) reproducibility of the technique, it is now possible to apply electroporation techniques to applications previously considered too sensitive for electroporation methodologies.
- The finer control over the delivered energy offers specific and important advantages for CULTURED EMBRYOS electroporation. As the thrust of NEPA21 protocols is to minimise delivered energy, this means that the targets are electroporated with less current (than competing device protocols).
- Only delivering the required energy (and no more) to porate the membrane is of utmost importance for viability post electroporation.
- The success of the NEPA21 for retina electroporation is evident by the Application and Publication information following.
- The NEPA21 system is supported by a suite of over 250 different electrode configurations, which further enhance the applicability of the system and empower researchers with further experimental options and opportunities.

APPLICATIONS


Transfection into Mouse-Rat: CULTURED EMBRYOS by Ex Utero Electroporation

Electroporation for mammalian embryos in the whole embryo culture system

Procedure:

1. Pre-culture embryos in the whole embryo culture system for 1.5-2 hours prior to electroporation.
2. Place the embryo in a Petri dish with Tyrode's solution.
3. Inject 0.1-0.5 µl of plasmid DNA into the brain ventricle with a fine capillary.
4. Apply square pulses using the electroporator (NEPA21/CUY21) and electrodes (chamber-type and forceps-type electrodes are available). 70V, 50 msec at 1-second interval, five pulses are applied for E10.5 mouse embryos.
5. Culture the electroporated embryos for 24-48 hours in the whole embryo culture system. Transfection of a fluorescent protein-expression vector into the developing rat cortex

Transfection of a fluorescent protein-expression vector into the developing rat cortex.




(A) EGFP-expression vector was electroporated into E11.5 rat telencephalon. The electroporated embryo was cultured in the whole embryo culture system (WEC).

(B) 24 hours after electroporation, EGFP-expression was specifically detected at the dorsal part of the telencephalon.

Tel: telencephalon;
Di: diencephalon

Transfection of fluorescent protein-expression vectors into the rat spinal cord.



PUBLICATIONS

Transfection into Mouse-Rat: CULTURED EMBRYOS by Ex Utero Electroporation

Ex utero mouse embryogenesis from pre-gastrulation to late organogenesis

Aguilera-Castrejon A, Oldak B, Shani T, Ghanem N, Itzkovich C, Slomovich S, Tarazi S, Bayerl J, Chugaeva V, Ayyash M, Ashoukhi S, Sheban D, Livnat N, Lasman L, Viukov S, Zerbib M, Addadi Y, Rais Y, Cheng S, Stelzer Y, Keren-Shaul H, Shlomo R, Massarwa R, Novershtern N, Maza I, Hanna JH.
Nature. 2021 May;593(7857):119-124.

IgSF11 homophilic adhesion proteins promote layer-specific synaptic assembly of the cortical interneuron subtype

Hayano Y, Ishino Y, Hyun JH, Orozco CG, Steinecke A, Potts E, Oisi Y, Thomas CI, Guerrero-Given D, Kim E, Kwon HB, Kamasawa N, Taniguchi H.
Sci Adv. 2021 Jul 14;7(29):eabf1600.

Autotaxin/lysophospholipase D-mediated lysophosphatidic acid signaling is required to form distinctive large lysosomes in the visceral endoderm cells of the mouse yolk sac.

Koike S, Keino-Masu K, Ohto T, Sugiyama F, Takahashi S, Masu M.
J Biol Chem. 2009 Nov 27;284(48):33561-70.

FGF8 signaling patterns the telencephalic midline by regulating putative key factors of midline development.

Okada T, Okumura Y, Motoyama J, Ogawa M.
Dev Biol. 2008 Aug 1;320(1):92-101.

Isolation and Characterization of Vasohibin-2 as a Homologue of VEGF-Inducible Endothelium-Derived Angiogenesis Inhibitor Vasohibin

Shibuya T, Watanabe K, Yamashita H, Shimizu K, Miyashita H, Abe M, Moriya T, Ohta H, Sonoda H, Shimosegawa T, Tabayashi K, Sato Y.
Arterioscler Thromb Vasc Biol. 2006 May;26(5):1051-7.

Misrouting of mitral cell progenitors in the Pax6/small eye rat telencephalon

Nomura T, Osumi N.

Development. 2004 Feb;131(4):787-96.

Manipulating gene expressions by electroporation in the developing brain of mammalian embryos

Takahashi M, Sato K, Nomura T, Osumi N.

Differentiation. 2002 Jun;70(4-5):155-62.

Pax6 regulates specification of ventral neurone subtypes in the hindbrain by establishing progenitor domains

Takahashi M, Osumi N.

Development . 2002 Mar;129(6):1327-38.

Gene Transfer into Cultured Mammalian Embryos by Electroporation

Osumi N, Inoue T.

Methods . 2001 May;24(1):35-42.