

Transfection into Mouse-Rat: **EMBRYOS IN UTERO** by Electroporation

The NEPA21 is the only device on the market to approach Electroporation from the perspective of optimising delivered energy.

- Compared to devices from other suppliers, the NEPA21 system offers the researcher a level of previously unavailable control over energy delivery to the electroporation target. This control is generated via unique electroporation pulse-output configurations, client-confirmed protocols and application-customised electrodes.
- With this market-leading control and (user-independent) reproducibility of the technique, it is now possible to apply electroporation techniques to applications previously considered too sensitive for electroporation methodologies.
- The finer control over the delivered energy offers specific and important advantages for EMBRYOS IN UTERO electroporation. As the thrust of NEPA21 protocols is to minimise delivered energy, this means that the targets are electroporated with less current (than competing device protocols).
- Only delivering the required energy (and no more) to porate the membrane is of utmost importance for viability post electroporation.
- The success of the NEPA21 for retina electroporation is evident by the Application and Publication information following.
- The NEPA21 system is supported by a suite of over 250 different electrode configurations, which further enhance the applicability of the system and empower researchers with further experimental options and opportunities.

APPLICATIONS

Transfection into Mouse-Rat: **EMBRYOS IN UTERO** by Electroporation

Gene transfer into embryonic brains using in utero electroporation technique


1) Materials



- In Vivo Electroporator: NEPA21/CUY21 SC/CUY21 EDIT (Nepa Gene Co., Ltd.)
- In Vivo Electrode: CUY650P3/CUY650P5 (Tweezers w/3mm or 5mm diameter platinum disk electrode, Nepa Gene Co., Ltd.)
- Aspirator tube assembly (Drummond)
- Optical fiber light (Technolight, Kenko, #KTS-100RSV)
- Sterile gauze (K-Pine, 7.5cm x 7.5cm)
- Surgical instruments: Fine forceps x 2, Surgery scissors x 2, Ring forceps, Needle holder, Surgical tape
- Nylon suture (Nesco, #HT1605NA75)
- Silk suture (D&G, #112451)



2) In Utero Electroporation Procedure




1. A 2 cm midline incision is then made in the abdominal wall along the linea alba using a set of forceps and scissors. A piece of sterile gauze with a hole cut in the center is placed over the incision, and one uterine horn is drawn out through the hole in the gauze.
2. After observing the orientation of the embryos through the wall of the uterine horn, a micropipette is inserted into the lateral ventricle, and 2-5 μ l of plasmid solution is injected by expiratory pressure using the aspirator tube assembly. When CAG-EGFP is used, a concentration of 1 μ g/ μ l is sufficient to visualize the migrating neurons.
3. After the injection, DNA solution containing 0.01% FastGreen can be seen through the uterine wall (red arrow).
4. After soaking the uterine horn with PBS, the head of embryo is pinched with a tweezers-type electrode, and electronic pulses are applied with the electroporator.

Electroporation settings for ICR mice

Age	Electrode (diameter)	Voltage	Pulse On	Pulse Off	Number of Pulses
E12.5	3mm	33V	30msec	970ms	4
E13.5	5mm	30V	50msec	950ms	4
E14.5	5mm	33V	50msec	950ms	4
E15-	5mm	35V	50msec	950ms	4

If viability was prioritized over transfection efficiency, the number of pulses can be changed to two. The actual current is displayed on the screen of the electroporator (NEPA21/CUY21/CUY21E, Nepa Gene). Make sure that the current would become 30-60mA. The current varies according to how the electrode applied or wetness of the uterine horn. Examine the gap between electrodes and the electrode contact areas to fit the current in the appropriate range. If the current is still above the range after the examination, change the voltage setting.

3) GFP expression



CAG-EGFP was injected into the both lateral ventricles of E14.5 mouse embryos and electronic pulses (33V, 50msec) were charged four times. 3 days later, the embryos (E17.5) were fixed and the brains were removed and examined under a fluorescence stereomicroscope (Fig. A).

Fluorescence was observed in the lateral region of the hemisphere onto which the anode had been placed and in the medial region of the opposite hemisphere.

Brains were frozen and sliced and the fluorescent image was obtained with a confocal laser microscope (Fig. B). GFP positive cells into which DNA was transferred at the ventricular zone (VZ) migrated to the intermediate zone (IZ) and cortical plate (CP). The arrowheads show the border between VZ and IZ and the border between IZ and CP. Dashed line show the border of tissues. VZ: Ventricular Zone, IZ: Intermediate Zone, CP: Cortical Plate

Hidegori Tabata and Kazunori Nakajima, Department of Anatomy, Keio University School of Medicine

PUBLICATIONS

Transfection into Mouse-Rat: **EMBRYOS IN UTERO** by Electroporation

E13.5 / E14.5

The microcephaly-associated transcriptional regulator AUTS2 cooperates with Polycomb complex PRC2 to produce upper-layer neurons in mice

Shimaoka K, Hori K, Miyashita S, Inoue YU, Tabe NKN, Sakamoto A, Hasegawa I, Nishitani K, Yamashiro K, Egusa SF, Tatsumoto S, Go Y, Abe M, Sakimura K, Inoue T, Imamura T, Hoshino M.

EMBO J. 2025 Jan 15.

Mouse / E14.0 / E15.5 / Lateral ventricle

Assembly of neuron- and radial glial-cell-derived extracellular matrix molecules promotes radial migration of developing cortical neurons

Mubuchi A, Takechi M, Nishio S, Matsuda T, Itoh Y, Sato C, Kitajima K, Kitagawa H, Miyata S.

Elife. 2024 Mar 21;12:RP92342.

Mouse / E12.5_E13.0 / E14.0_E15.0 / E16.0_E17.0 / Cerebral cortex / Parietal lobe Somatosensory area**Analyses of Conditional Knockout Mice for Pogz, a Gene Responsible for Neurodevelopmental Disorders in Excitatory and Inhibitory Neurons in the Brain**

Hamada N, Nishijo T, Iwamoto I, Shifman S, Nagata KI.

Cells. 2024 Mar 19;13(6):540.

Mouse / Lateral_ventricle / E14.5**Lmo4 synergizes with Fezf2 to promote direct in vivo reprogramming of upper layer cortical neurons and cortical glia towards deep-layer neuron identities**

Felske T, Tocco C, Péron S, Harb K, Alfano C, Galante C, Berninger B, Studer M.

PLoS Biol. 2023 Aug 8;21(8):e3002237.

Mouse / E14.0 / Lateral ventricle**Histological Analysis of a Mouse Model of the 22q11.2 Microdeletion Syndrome**

Tabata H, Mori D, Matsuki T, Yoshizaki K, Asai M, Nakayama A, Ozaki N, Nagata KI.

Biomolecules. 2023 Apr 27;13(5):763.

Mouse / E13.5 / E14.5 / Cerebral ventricles**p53/p21 pathway activation contributes to the ependymal fate decision downstream of GemC1**

Ortiz-Álvarez G, Fortoul A, Srivastava A, Moreau MX, Bouloudi B, Mailhes-Hamon C, Delgehyr N, Faucourt M, Bahin M, Blugeon C, Breau M, Géli V, Causeret F, Meunier A, Spassky N.

Cell Rep. 2022 Dec 13;41(11):111810.

Mouse / Cerebral cortex / Parietal lobe Somatosensory area / E14.0**Variant-specific changes in RAC3 function disrupt corticogenesis in neurodevelopmental phenotypes**

Scala M, Nishikawa M, Ito H, Tabata H, Khan T, Accogli A, Davids L, Ruiz A, Chiurazzi P, Cericola G, Schulte B, Monaghan KG, Begtrup A, Torella A, Pinelli M, Denommé-Pichon AS, Vitobello A, Racine C, Mancardi MM, Kiss C, Guerin A, Wu W, Gabau Vila E, Mak BC, Martinez-Agosto JA, Gorin MB, Duz B, Bayram Y, Carvalho CMB, Vengoechea JE, Chitayat D, Tan TY, Callewaert B, Kruse B, Bird LM, Faivre L, Zollino M, Biskup S; Undiagnosed Diseases Network; Telethon Undiagnosed Diseases Program; Striano P, Nigro V, Severino M, Capra V, Costain G, Nagata KI.

Brain. 2022 Sep 14;145(9):3308-3327.

Mouse / Cerebral cortex / E14.5**Tbr1 Misexpression Alters Neuronal Development in the Cerebral Cortex**

Crespo I, Pignatelli J, Kinare V, Méndez-Gómez HR, Esgleas M, Román MJ, Canals JM, Tole S, Vicario C.

Mol Neurobiol. 2022 Sep;59(9):5750-5765.

Mouse**Heterogeneous nuclear ribonucleoprotein U (HNRNPU) safeguards the developing mouse cortex**

Sapir T, Kshirsagar A, Gorelik A, Olender T, Porat Z, Scheffer IE, Goldstein DB, Devinsky O, Reiner O.

Nat Commun. 2022 Jul 21;13(1):4209.

Mouse / E13.5 / E14.5 / Lateral ventricle**Endosomal trafficking defects alter neural progenitor proliferation and cause microcephaly**

Carpentieri JA, Di Cicco A, Lampic M, Andreau D, Del Maestro L, El Marjou F, Coquand L, Bahi-Buisson N, Brault JB, Baffet AD.

Nat Commun. 2022 Jan 10;13(1):16.

Mouse / E14.0 / Lateral ventricle**Heterogeneous nuclear ribonucleoprotein U (HNRNPU) safeguards the developing mouse cortex**

Sapir T, Kshirsagar A, Gorelik A, Olender T, Porat Z, Scheffer IE, Goldstein DB, Devinsky O, Reiner O.
Nat Commun. 2022 Jul 21;13(1):4209.

Mouse / E14.0 / Cerebral cortex / Parietal lobe Somatosensory area**Impaired Function of PLEKHG2, a Rho-Guanine Nucleotide-Exchange Factor, Disrupts Corticogenesis in Neurodevelopmental Phenotypes**

Nishikawa M, Ito H, Tabata H, Ueda H, Nagata KI.
Cells. 2022 Feb 16;11(4):696.

Mouse / Lateral ventricle**Simultaneous two-photon imaging of action potentials and subthreshold inputs in vivo**

Bando Y, Wenzel M, Yuste R.
Nat Commun. 2021 Dec 10;12(1):7229

Mouse / E15.0 / Cerebral cortex / Neuroprogenitor cells**IgSF11 homophilic adhesion proteins promote layer-specific synaptic assembly of the cortical interneuron subtype**

Hayano Y, Ishino Y, Hyun JH, Orozco CG, Steinecke A, Potts E, Oisi Y, Thomas CI, Guerrero-Given D, Kim E, Kwon HB, Kamasawa N, Taniguchi H.
Sci Adv. 2021 Jul 14;7(29):eabf1600.

Mouse / E14.5 / Cerebral cortex**A pericellular hyaluronan matrix is required for the morphological maturation of cortical neurons**

Takechi M, Oshima K, Nadano D, Kitagawa H, Matsuda T, Miyata S.
Biochim Biophys Acta Gen Subj. 2020 Oct;1864(10):129679.

Mouse / E13.5 / Cerebral cortex**Huntington's disease alters human neurodevelopment**

Barnat M, Capizzi M, Aparicio E, Boluda S, Wennagel D, Kacher R, Kassem R, Lenoir S, Agasse F, Braz BY, Liu JP, Ighil J, Tessier A, Zeitlin SO, Duyckaerts C, Dommergues M, Durr A, Humbert S.
Science. 2020 Aug 14;369(6505):787-793.

Mouse / Neural tube**In utero gene transfer system for embryos before neural tube closure reveals a role for Hmga2 in the onset of neurogenesis**

Naohiro Kuwayama, Yusuke Kishi, Yurie Maeda, Yurie Nishiumi, Yutaka Suzuki, Haruhiko Koseki, Yusuke Hirabayashi, Yukiko Gotoh
bioRxiv May 15, 2020

Mouse / E14.5 / Cerebral cortex / Neuroprogenitor cells**Knockdown of Son, a mouse homologue of the ZTTK syndrome gene, causes neuronal migration defects and dendritic spine abnormalities**

Ueda M, Matsuki T, Fukada M, Eda S, Toya A, Iio A, Tabata H, Nakayama A.
Mol Brain. 2020 May 24;13(1):80.

Mouse / E15.5 / Lateral ventricle**DNA Repair Protein RAD51 Enhances the CRISPR/Cas9-mediated Knock-In Efficiency in Brain Neurons**

Taiga Kurihara, Emi Kouyama-Suzuki, Michiru Satoga, Xue Li, Moataz Badawi, Thiha, Deeba Noreen Baig, Toru Yanagawa, Takeshi Uemura, Takuma Mori, Katsuhiko Tabuchi

Biochem Biophys Res Commun. 2020 Apr 9;524(3):621-628.

Mouse / E14.5_E15.5 / Lateral ventricle**In Vivo Imaging of the Coupling Between Neuronal and CREB Activity in the Mouse Brain**

Tal Laviv, Benjamin Scholl, Paula Parra-Bueno, Beth Foote, Chuqiu Zhang, Long Yan, Yuki Hayano, Jun Chu, Ryohei Yasuda

Neuron, 105 (5), 799-812.e5 2020 Mar 4

Mouse / Otocyst**Prenatal Electroporation-Mediated Gene Transfer Restores Slc26a4 Knock-Out Mouse Hearing and Vestibular Function**

Hiroki Takeda, Toru Miwa, Min Young Kim, Byung Yoon Choi, Yori-hisa Orita, Ryosei Minoda

Sci Rep, 9 (1), 17979 2019 Nov 29

Mouse / E15.0 / Postnatal electroporation / P2 / Adult Mouse electroporation / P21**Primary cilium-dependent cAMP/PKA signaling at the centrosome regulates neuronal migration**

Stoufflet J, Chaulet M, Doulazmi M, Fouquet C, Dubacq C, Métin C, Schneider-Maunoury S, Trembleau A, Vincent P, Caillé I.

Sci Adv. 2020 Sep 2;6(36):eaba3992

Mouse / E14.5 / Lateral ventricle**A Novel LGI1 Missense Mutation Causes Dysfunction in Cortical Neuronal Migration and Seizures**

Feng Liu, Chao Du, Xin Tian, Yuanlin Ma, Bei Zhao, Yin Yan, Zijun Lin, Peijia Lin, Ruijiao Zhou, Xuefeng Wang

Brain Res, 1721, 146332 2019 Oct 15

Mouse / Embryonic inner ear / Otocyst**Role of Dach1 Revealed Using a Novel Inner Ear-Specific Dach1-knockdown Mouse Model**

Toru Miwa, Ryosei Minoda, Yoshihide Ishikawa, Tomohito Kajii, Yori-hisa Orita, Takahiro Ohyama

Biol Open, 8 (8) 2019 Aug 20

Mouse / E14.5 / Lateral ventricle**Estrogen Receptor α Promotes Cav1.2 Ubiquitination and Degradation in Neuronal Cells and in APP/PS1 Mice**

Yu-Jie Lai, Bing-Lin Zhu, Fei Sun, Dong Luo, Yuan-Lin Ma, Bio Luo, Jing Tang, Ming-Jian Xiong, Lu Liu, Yan Long, Xiao-Tong Hu, Ling He, Xiao-Juan Deng, John H Zhang, Jian Yang, Zhen Yan, Guo-Jun Chen

Aging Cell, 18 (4), e12961 Aug 2019

Mouse In Utero / E14.5 / Cerebral ventricles**In Vivo Single-Cell Genotyping of Mouse Cortical Neurons Transfected With CRISPR/Cas9**

Steinecke A, Kurabayashi N, Hayano Y, Ishino Y, Taniguchi H.

Cell Rep, 28 (2), 325-331.e4 2019 Jul 9

Mouse / E14.5 / Lateral ventricle**Both Excitatory and Inhibitory Neurons Transiently Form Clusters at the Outermost Region of the Developing Mammalian Cerebral Neocortex**

Shin M, Kitazawa A, Yoshinaga S, Hayashi K, Hirata Y, Dehay C, Kubo KI, Nakajima K.
J Comp Neurol. 2019 Jul 1;527(10):1577-1597.

Mouse / E13.5 / Cerebral cortex / Somatosensory cortex**Pathological mTOR Mutations Impact Cortical Development**

Tarkowski B, Kuchcinska K, Blazejczyk M, Jaworski J.
Hum Mol Genet. 2019 Jul 1;28(13):2107-2119.

Mouse / E15.0 / Cerebral cortex / Glucocorticoid receptor**Persistence of Learning-Induced Synapses Depends on Neurotrophic Priming of Glucocorticoid Receptors**

Arango-Lievano M, Borie AM, Dromard Y, Murat M, Desarmenien MG, Garabedian MJ, Jeanneteau F.
Proc Natl Acad Sci U S A. 2019 Jun 25;116(26):13097-13106.

Mouse / Cerebral cortex /**Role of Per3, a Circadian Clock Gene, in Embryonic Development of Mouse Cerebral Cortex**

Mariko Noda, Ikuko Iwamoto, Hidenori Tabata, Takanori Yamagata, Hidenori Ito, Koh-Ichi Nagata Sci Rep, 9 (1), 5874 2019 Apr 10.

Mouse / E14.5 / Lateral ventricles**MicroRNA-129-5p Is Regulated by Choline Availability and Controls EGF Receptor Synthesis and Neurogenesis in the Cerebral Cortex**

Trujillo-Gonzalez I, Wang Y, Friday WB, Vickers KC, Toth CL, Molina-Torres L, Surzenko N, Zeisel SH.
FASEB J. 2019 Mar;33(3):3601-3612.

Mouse / E12.5 / E15.5 / E16.0 / Cortical progenitors**Intersectional Monosynaptic Tracing for Dissecting Subtype-Specific Organization of GABAergic Interneuron Inputs**

Yetman MJ, Washburn E, Hyun JH, Osakada F, Hayano Y, Zeng H, Callaway EM, Kwon HB, Taniguchi H.
Nat Neurosci. 2019 Mar;22(3):492-502.

Single cell profiling of CRISPR/Cas9-induced OTX2 deficient retinas reveals fate switch from restricted progenitors

Miruna G. Ghinia Tegla, Diego F. Buenaventura, Diana Y. Kim, Cassandra Thakurdin, Kevin C. Gonzalez, M. Emerson
bioRxiv February 02, 2019

An Optogenetic Approach to Studies of the Mechanisms of Heterosynaptic Plasticity in Neocortical Neurons

N. A. Simonova, N. V. Bal, P. M. Balaban, M. A. Volgushev, A. Y. Malyshev
Neuroscience and Behavioral Physiology volume 49, pages208–215(2019)

Drebrin-like (Dbnl) Controls Neuronal Migration via Regulating N-Cadherin Expression in the Developing Cerebral Cortex

Seika Inoue, Kanehiro Hayashi, Kyota Fujita, Kazuhiko Tagawa, Hitoshi Okazawa, Ken-Ichiro Kubo, Kazunori Nakajima
J Neurosci, 39 (4), 678-691 2019 Jan 23

Comparative Evaluation of Genetically Encoded Voltage Indicators

Yuki Bando, Masayuki Sakamoto, Samuel Kim, Inbal Ayzenshtat, Rafael Yuste
Cell Rep, 26 (3), 802-813.e4 2019 Jan 15

Niche-derived laminin-511 Promotes Midbrain Dopaminergic Neuron Survival and Differentiation Through YAP

Dawei Zhang, Shanzheng Yang, Enrique M Toledo, Daniel Gyllborg, Carmen Saltó, J Carlos Villaescusa, Ernest Arenas

Sci Signal, 10 (493) 2017 Aug 22

PLGF, a placental marker of fetal brain defects after in utero alcohol exposure.

Lecuyer M, Laquerrière A, Bekri S, Lesueur C, Ramdani Y, Jégou S, Uguen A, Marcorelles P, Marret S, Gonzalez BJ

Acta Neuropathol Commun. 2017 Jun 6;5(1):44.

Developmental activities of the complement pathway in migrating neurons.

Gorelik A1, Sapir T, Haffner-Krausz R, Olender T, Woodruff TM, Reiner O

Nat Commun. 2017 May 2;8:15096.

PLGF, a placental marker of fetal brain defects after in utero alcohol exposure.

Lecuyer M, Laquerrière A, Bekri S, Lesueur C, Ramdani Y, Jégou S, Uguen A, Marcorelles P, Marret S, Gonzalez BJ

Acta Neuropathol Commun. 2017 Jun 6;5(1):44.

In vivo genome editing via CRISPR/Cas9 mediated homology-independent targeted integration.

Suzuki K, Tsunekawa Y, Hernandez-Benitez R, Wu J, Zhu J, Kim EJ, Hatanaka F, Yamamoto M, Araoka T, Li Z, Kurita M, Hishida T, Li M, Aizawa E, Guo S, Chen S, Goebel A, Soligalla RD, Qu J, Jiang T, Fu X, Jafari M, Esteban CR, Berggren WT, Lajara J, Nuñez-Delicado E, Guillen P, Campistol JM, Matsuzaki F, Liu GH, Magistretti P, Zhang K, Callaway EM, Zhang K, Belmonte JC

Nature. 2016 Dec 1;540(7631):144-149

Early postnatal GABAA receptor modulation reverses deficits in neuronal maturation in a conditional neurodevelopmental mouse model of DISC1.

Saito A, Taniguchi Y, Rannals MD, Merfeld EB, Ballinger MD, Koga M, Ohtani Y, Gurley DA, Sedlak TW, Cross A, Moss SJ, Brandon NJ, Maher BJ, Kamiya A

Mol Psychiatry. 2016 Oct;21(10):1449-59.

CRISPR/Cas9-mediated Gene Knock-Down in Post-Mitotic Neurons

Christoph Straub, Adam J Granger, Jessica L Saulnier, Bernardo L Sabatini

PLoS One, 9 (8), e105584 2014 Aug 20 eCollection 2014

Targeted DNA demethylation in vivo using dCas9-peptide repeat and scFv-TET1 catalytic domain fusions.

Morita S, Noguchi H, Horii T, Nakabayashi K, Kimura M, Okamura K, Sakai A, Nakashima H, Hata K, Nakashima K, Hatada I

Nat Biotechnol. 2016 Aug 29 /nbt.3658. [Epub ahead of print]

High-Throughput, High-Resolution Mapping of Protein Localization in Mammalian Brain by In Vivo Genome Editing.

Mikuni T, Nishiyama J, Sun Y, Kamasawa N, Yasuda R

Cell. 2016 Jun 16;165(7):1803-17.

Targeted in vivo genetic manipulation of the mouse or rat brain by in utero electroporation with a triple-electrode probe.

Szczurkowska J, Cwetsch AW, dal Maschio M, Ghezzi D, Ratto GM, Cancedda L

Nat Protoc. 2016 Mar;11(3):399-412.

Developing a De Novo Targeted Knock-In Method Based on in Utero Electroporation Into the Mammalian Brain

Yuji Tsunekawa, Raymond Kunikane Terhune, Ikumi Fujita, Atsunori Shitamukai, Taeko Suetsugu, Fumio Matsuzaki

Development, 143 (17), 3216-22 2016 Sep 1

Radmis, a Novel Mitotic Spindle Protein that Functions in Cell Division of Neural Progenitors.

Yumoto T, Nakadate K, Nakamura Y, Sugitani Y, Sugitani-Yoshida R, Ueda S, Sakakibara S

PLoS One. 2013 Nov 8;8(11):e79895.

Screening for candidate genes involved in the production of mouse subventricular zone proliferative cells and an estimation of their changes in evolutionary pressure during primate evolution.

Tabata H, Hachiya T, Nagata K, Sakakibara Y, Nakajima K.

Front Neuroanat. 2013 Jul 31;7:24.

Development of cerebellar neurons and glias revealed by in utero electroporation: Golgi-like labeling of cerebellar neurons and glias.

Kita Y, Kawakami K, Takahashi Y, Murakami F.

PLoS One. 2013 Jul 23;8(7):e70091.

Shootin1 Acts in Concert with KIF20B to Promote Polarization of Migrating Neurons.

Sapir T, Levy T, Sakakibara A, Rabinkov A, Miyata T, Reiner O.

J Neurosci. 2013 Jul 17;33(29):11932-48

Selective and regulated gene expression in murine Purkinje cells by in utero electroporation.

Nishiyama J, Hayashi Y, Nomura T, Miura E, Kakegawa W, Yuzaki M.

Eur J Neurosci. 2012 Oct;36(7):2867-76.

Persistent cortical plasticity by upregulation of chondroitin 6-sulfation.

Miyata S, Komatsu Y, Yoshimura Y, Taya C, Kitagawa H.

Nat Neurosci. 2012 Jan 15;15(3):414-22, S1-2.

Imaging brain electric signals with genetically targeted voltage-sensitive fluorescent proteins

Akemann W, Mutoh H, Perron A, Rossier J, Knöpfel T.

Nature Methods, Volume 7, Number8 Pages 643-649, August 2010

Protein phosphatase 2A facilitates axonogenesis by dephosphorylating CRMP2.

Zhu LQ, Zheng HY, Peng CX, Liu D, Li HL, Wang Q, Wang JZ.

J Neurosci. 2010 Mar 10;30(10):3839-48.

Knockdown of DISC1 by in utero gene transfer disturbs postnatal dopaminergic maturation in the frontal cortex and leads to adult behavioral deficits.

Niwa M, Kamiya A, Murai R, Kubo K, Gruber AJ, Tomita K, Lu L, Tomisato S, Jaaro-Peled H, Seshadri S, Hiyama H, Huang B, Kohda K, Noda Y, O'Donnell P, Nakajima K, Sawa A, Nabeshima T.

Neuron. 2010 Feb 25;65(4):480-9.

Neurotrophin-3 is involved in the formation of apical dendritic bundles in cortical layer 2 of the rat.

Miyashita T, Wintzer M, Kurotani T, Konishi T, Ichinohe N, Rockland KS.

Cereb Cortex. 2010 Jan;20(1):229-40.

Involvement of the cytoplasmic C-terminal domain of connexin43 in neuronal migration.

Cina C, Maass K, Theis M, Willecke K, Bechberger JF, Naus CC.

J Neurosci. 2009 Feb 18;29(7):2009-21.

Electroporation-mediated gene transfer to the developing mouse inner ear.

Brigande JV, Gubbels SP, Woessner DW, Jungwirth JJ, Bresee CS.
Methods Mol Biol. 2009;493:125-39.

Labeling embryonic mouse central nervous system cells by in utero electroporation.

Tabata H, Nakajima K
Dev Growth Differ. 2008 Aug;50(6):507-11.

Sirt1 contributes critically to the redox-dependent fate of neural progenitors

Prozorovski T, Schulze-Topphoff U, Glumm R, Baumgart J, Schröter F, Ninnemann O, Siegert E, Bendix I, Brüstle O, Nitsch R, Zipp F, Aktas O.
Nature Cell Biology, Volume 10, Number 4, Pages 385-394, April 2008

RECK modulates Notch signaling during cortical neurogenesis by regulating ADAM10 activity

Muraguchi T, Takegami Y, Ohtsuka T, Kitajima S, Chandana EP, Omura A, Miki T, Takahashi R, Matsumoto N, Ludwig A, Noda M, Takahashi C.
Nature Neuroscience, Volume 10, Number 7, Pages 838-845, July 2007

Direct evidence that ventral forebrain cells migrate to the cortex and contribute to the generation of cortical myelinating oligodendrocytes.

Nakahira E, Kagawa T, Shimizu T, Goulding MD, Ikenaka K
Dev Biol. 2006 Mar 1;291(1):123-31

Radial Migration of Superficial Layer Cortical Neurons Controlled by Novel Ig Cell Adhesion Molecule MDGA1

Takeuchi A, O'Leary DD.
The Journal of Neuroscience, Volume 26, Issue 17, Pages 4460-4464, 26 April 2006

In vivo electroporation in the embryonic mouse central nervous system

Saito T
Nature Protocols, Volume 1, Number 3, Pages 1552-1558, 2006

G Protein betagamma Subunits and AGS3 Control Spindle Orientation and Asymmetric Cell Fate of Cerebral Cortical Progenitors

Sanada K, Tsai LH.
Cell. 2005 Jul 15;122(1):119-31.

Birth-date-dependent segregation of the mouse cerebral cortical neurons in reaggregation cultures

Ajioka I, Nakajima K.
Eur J Neurosci. 2005 Jul;22(2):331-42.

Commissural neuron identity is specified by a homeodomain protein, Mbh1, that is directly downstream of Math1

Saba R, Johnson JE, Saito T.
Development. 2005 May;132(9):2147-55.

Neuronal generation, migration, and differentiation in the mouse hippocampal primordium as revealed by enhanced green fluorescent protein gene transfer by means of in utero electroporation

Nakahira E, Yuasa S.
J Comp Neurol. 2005 Mar 14;483(3):329-40.

A NUDEL-dependent mechanism of neurofilament assembly regulates the integrity of CNS neurons

Nguyen MD, Shu T, Sanada K, Larivière RC, Tseng HC, Park SK, Julien JP, Tsai LH.
Nat Cell Biol. 2004 Jul;6(7):595-608.

Hes genes regulate size, shape and histogenesis of the nervous system by control of the timing of neural stem cell differentiation

Hatakeyama et al.

Development, Volume 131, Issue 22, Pages 5539-5550, November 2004

Mouse / E11.5**Generation of Reelin-Positive Marginal Zone Cells from the Caudomedial Wall of Telencephalic Vesicles**

Takiguchi-Hayashi K, Sekiguchi M, Ashigaki S, Takamatsu M, Hasegawa H, Suzuki-Migishima R, Yokoyama M, Nakanishi S, Tanabe Y.

J Neurosci. 2004 Mar 3;24(9):2286-95.

Mouse / E13.0**Multipolar Migration: The Third Mode of Radial Neuronal Migration in the Developing Cerebral Cortex**

Tabata H, Nakajima K.

J Neurosci. 2003 Nov 5;23(31):9996-10001.

Mouse / E14.0**The in vivo roles of STEF/Tiam1, Rac1 and JNK in cortical neuronal migration**

Kawauchi T, Chihama K, Nabeshima Y, Hoshino M.

EMBO J. 2003 Aug 15;22(16):4190-201.

Layering defect in p35 deficiency is linked to improper neuronal-glia interaction in radial migration

Gupta et al.

Nature Neuroscience, Volume 6, Number 12, Pages 1284-1291, December 2003

Mouse / E14.0**Serine 732 Phosphorylation of FAK by Cdk5 Is Important for Microtubule Organization, Nuclear Movement, and Neuronal Migration**

Xie Z, Sanada K, Samuels BA, Shih H, Tsai LH.

Cell. 2003 Aug 22;114(4):469-82.

Mouse / E11.5 / Spinal cord**Mammalian BarH1 Confers Commissural Neuron Identity on Dorsal Cells in the Spinal Cord**

Saba R, Nakatsuji N, Saito T.

J Neurosci. 2003 Mar 15;23(6):1987-91.

Emx2 patterns the neocortex by regulating FGF positional signaling

Fukuchi-Shimogori T, Grove EA.

Nat Neurosci. 2003 Aug;6(8):825-31.

Roles of the Basic Helix-Loop-Helix Genes Hes1 and Hes5 in Expansion of Neural Stem Cells of the Developing Brain

Ohtsuka et al.

The Journal of Biological Chemistry, Volume 276, Issue 32, Pages 30467-30474, 10 August 2001

Mouse / E12.5_E13 / E14_E15 / E16_E17**Efficient in utero gene transfer system to the developing mouse brain using electroporation: visualization of neuronal migration in the developing cortex**

Tabata H, Nakajima K.

Neuroscience. 2001;103(4):865-72.