

In utero electroporation parameter

Stage	Diameter	Electrode	Voltage	Pulse length	Pulse interval	Cycles	Current
E12.5	3mm	CUY650P3	33V	30 msec.	970 msec.	4	30-60mA
E13.5	5mm	CUY650P5	30V	50 msec.	950 msec.	4	30-60mA
E14.5	5mm	CUY650P5	33V	50 msec.	950 msec.	4	30-60mA
E15	5mm	CUY650P5	35V	50 msec.	950 msec.	4	30-60mA